




BID Antibody, FITC conjugated

Product Code	CSB-PA002698HC01HU
Uniprot NO.	P55957
Target Names	BID
Alternative Names	Apoptotic death agonist antibody; Apoptotic death agonist BID antibody; BH3 interacting domain death agonist antibody; BH3 interacting domain death agonist p11 antibody; BH3 interacting domain death agonist p13 antibody; BH3 interacting domain death agonist p15 antibody; BH3-interacting domain death agonist p11 antibody; BID antibody; BID isoform ES(1b) antibody; BID isoform L(2) antibody; BID isoform Si6 antibody; BID_HUMAN antibody; Desmocollin type 4 antibody; FP497 antibody; Human BID coding sequence antibody; MGC15319 antibody; MGC42355 antibody; p11 BID antibody; p13 BID antibody; p15 BID antibody; p22 BID antibody
Raised in	Rabbit
Species Reactivity	Human
Immunogen	Recombinant Human BH3-interacting domain death agonist protein (1-195AA)
Immunogen Species	Homo sapiens (Human)
Conjugate	FITC
Isotype	IgG
Purification Method	>95%, Protein G purified
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Image	 BID Antibody, FITC conjugated