





SMAD3 (Ab-425) Antibody

Product Code	CSB-PA794085
Uniprot NO.	P84022
Target Names	SMAD3
Alternative Names	DKFZP586N0721 antibody; DKFZp686J10186 antibody; hMAD 3 antibody; hMAD-3 antibody; hSMAD3 antibody; HSPC193 antibody; HST17436 antibody; JV15 2 antibody; JV15-2 antibody; JV152 antibody; LDS1C antibody; LDS3 antibody; MAD (mothers against decapentaplegic Drosophila) homolog 3 antibody; MAD homolog 3 antibody; Mad homolog JV15 2 antibody; Mad protein homolog antibody; MAD; mothers against decapentaplegic homolog 3 antibody; Mad3 antibody; MADH 3 antibody; MADH3 antibody; MGC60396 antibody; Mothers against decapentaplegic homolog 3 antibody; Mothers against DPP homolog 3 antibody; SMA and MAD related protein 3 antibody; SMAD 3 antibody; SMAD antibody; SMAD family member 3 antibody; SMAD; mothers against DPP homolog 3 antibody; Smad3 antibody; SMAD3_HUMAN antibody
Raised in	Rabbit
Species Reactivity	Human,Mouse,Rat
Immunogen	Peptide sequence around aa.423~427 (C-S-S-V-S) derived from Human Smad3.
Immunogen Species	Homo sapiens (Human)
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Concentration	It differs from different batches. Please contact us to confirm it.
Tested Applications	ELISA,IHC,IF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Image



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Smad3(Ab-425) Antibody(left) or the same antibody preincubated with blocking peptide(right).



Immunofluorescence staining of methanol-fixed HeLa cells using Smad3(Ab-425) Antibody.
