




Phospho-HTRA2 (S212) Antibody

Product Code	CSB-PA008882
Uniprot NO.	O43464
Target Names	HTRA2
Alternative Names	HTRA2; OMI; PRSS25; Serine protease HTRA2, mitochondrial; High temperature requirement protein A2; HtrA2; Omi stress-regulated endoprotease; Serine protease 25; Serine proteinase OMI
Raised in	Rabbit
Species Reactivity	Human,Mouse,Rat
Immunogen	Synthesized peptide derived from Human HtrA2 around the phosphorylation site of S212.
Immunogen Species	Homo sapiens (Human)
Conjugate	Non-conjugated
Isotype	IgG
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Tested Applications	IHC, ELISA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Image	 Phospho-HTRA2 (S212) Antibody