





ANGPT1 Antibody

Product Code	CSB-PA235182
Uniprot NO.	Q15389
Target Names	ANGPT1
Alternative Names	AGP 1 antibody; AGP1 antibody; AGPT antibody; ANG 1 antibody; ANG antibody; ANG-1 antibody; ANG1 antibody; Angiopoietin-1 antibody; Angiopoietin1 antibody; ANGPT1_HUMAN antibody; ANGPT 1 antibody; Angpt1 antibody; KIAA0003 antibody; OTTHUMP00000227592 antibody; OTTHUMP00000227599 antibody; OTTHUMP00000227600 antibody
Raised in	Rabbit
Species Reactivity	Human,Mouse,Rat
Immunogen	Fusion protein of Human ANGPT1
Immunogen Species	Homo sapiens (Human)
Conjugate	Non-conjugated
Isotype	IgG
Purification Method	Antigen affinity purification
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Tested Applications	ELISA,IHC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Image



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA235182(ANGPT1 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA235182(ANGPT1 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)