




PIK3R1 Antibody, HRP conjugated

Product Code	CSB-PA018003LB01HU
Uniprot NO.	P27986
Target Names	PIK3R1
Alternative Names	GRB1 antibody; p85 alpha antibody; p85 antibody; P85A_HUMAN antibody; Phosphatidylinositol 3 kinase associated p 85 alpha antibody; Phosphatidylinositol 3 kinase regulatory 1 antibody; Phosphatidylinositol 3 kinase; regulatory subunit; polypeptide 1 (p85 alpha) antibody; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alpha antibody; Phosphatidylinositol 3-kinase regulatory subunit alpha antibody; Phosphoinositide 3 kinase; regulatory subunit 1 (alpha) antibody; PI3 kinase p85 antibody; PI3 kinase p85 subunit alpha antibody; PI3-kinase regulatory subunit alpha antibody; PI3-kinase subunit p85-alpha antibody; PI3K antibody; PI3K p85 antibody; PI3K regulatory subunit alpha antibody; Pik3r1 antibody; PtdIns 3 kinase p85 alpha antibody; PtdIns-3-kinase regulatory subunit alpha antibody; PtdIns-3-kinase regulatory subunit p85-alpha antibody
Raised in	Rabbit
Species Reactivity	Human
Immunogen	Recombinant Human Phosphatidylinositol 3-kinase regulatory subunit alpha protein (160-376AA)
Immunogen Species	Homo sapiens (Human)
Conjugate	HRP
Isotype	IgG
Purification Method	>95%, Protein G purified
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Tested Applications	ELISA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Image	 PIK3R1 Antibody, HRP conjugated