




OR1D2 Antibody, FITC conjugated

Product Code	CSB-PA016490LC01HU
Uniprot NO.	P34982
Target Names	OR1D2
Alternative Names	OR1D2; OLFR1; Olfactory receptor 1D2; Olfactory receptor 17-4; OR17-4; Olfactory receptor OR17-6; Olfactory receptor-like protein HGMP07E
Raised in	Rabbit
Species Reactivity	Human
Immunogen	Recombinant Human Olfactory receptor 1D2 protein (159-196AA)
Immunogen Species	Homo sapiens (Human)
Conjugate	FITC
Isotype	IgG
Purification Method	>95%, Protein G purified
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Image	 OR1D2 Antibody, FITC conjugated