




Recombinant Human Sphingosine 1-phosphate receptor 3 (S1PR3)

Product Code	CSB-CF859932HU
Target Names	S1PR3
Uniprot NO.	Q99500
Species	Homo sapiens (Human)
Source	in vitro E.coli expression system
Expression Region	1-378
Target Protein Sequence	MATALPPRLQPVRGNETLREHYQYVGKLAGRLKEASEGSTLTTVLFLVICSFI VLENLMV LIAIWKNNKFHNRMYFFIGNLALCDLLAGIAYKVNILMSGKKTFSLSPTVWFL REGSMFV ALGASTCSLLAIAIERHLMIKMRPYDANKRHRVFLIGMCWLIAFTLGALPIL GWNCLH NLPDCSTILPLYSKKYIAFCISIFTAILVTIVILYARIYFLVKSSSRKVANHNNSE RSMA LLRTVVIVVSVFIACWSPLFILFLIDVACRVQACPILFKAQWFIVLAVLNSAMN PVIYTL ASKEMRRAFFRLVCNCLVRGRGARASPIQPALDPSRSKSSSSNNSSHSPKVKE DLPHTAP SSCIMDKNAALQNGIFCN
Protein Length	full length protein
Tag Info	N-terminal 10xHis-tagged
Image	



Recombinant Human Sphingosine 1-phosphate receptor 3 (S1PR3)