



Recombinant Sheep Thyrotropin receptor (TSHR), partial

Product Code	025131SH1
Purity	>85% (SDS-PAGE)
Target Names	TSHR
Uniprot NO.	P56495
Species	Ovis aries (Sheep)
Source	Yeast
Expression Region	22-413aa
Target Protein Sequence	ERCPSPPCECRQEDDFRVTCKDIQRIPLPSTQTLKFIETHLKTIPSRAFSNLPN ISRI YLSIDATLQQLESHSFYNLSKVTHIEIRNTRSLTYIDSGALKELPLLKFLGIFNT GLRVF PDLTKIYSTDVFFILEITDNPYMTSVPANAFQGLSNETLTLKLYNNGFTSIQGH AFNGTK LDAVYLNKNKYLTVIDQDAFAGVYSGPTLLDISYTSVTALPSKGLEHLKELIA RNTWTLK KLPLSLSFLHLTRADLSYP SHCCAFKNQKNIRGILQSLMCNESSIWGLRQRKS ASALNGP FYQEYEEDLGDGSAGYKENS KFQDTHSN SHYYVFFEDQEDEIIGFGQELKNP QEETLQAF DNHYDYTVCGGSEEMVCTPKSDEFNPCEDIMG
Protein Length	Partial
Tag Info	N-terminal His-tagged/Tag-Free
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.