




KIR2DL5A/KIR2DL5B Antibody

Product Code	CSB-PA944493
Uniprot NO.	Q8N109/Q8NHK3
Target Names	KIR2DL5A/KIR2DL5B
Raised in	Rabbit
Species Reactivity	Human
Immunogen	Synthetic peptide from Human protein at AA range: 121-170
Immunogen Species	Homo sapiens (Human)
Conjugate	Non-conjugated
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Tested Applications	IHC,ELISA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Image	



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200