



GPR85 Antibody, HRP conjugated

Product Code	CSB-PA009837LB01HU
Uniprot NO.	P60893
Target Names	GPR85
Alternative Names	GPR85; SREB2; Probable G-protein coupled receptor 85; Super conserved receptor expressed in brain 2
Raised in	Rabbit
Species Reactivity	Human
Immunogen	Recombinant Human Probable G-protein coupled receptor 85 protein (229-267AA)
Immunogen Species	Homo sapiens (Human)
Conjugate	HRP
Isotype	IgG
Purification Method	>95%, Protein G purified
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Tested Applications	ELISA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Image	

GPR85 Antibody, HRP conjugated

