




UBQLN2 Antibody, FITC conjugated

Product Code	CSB-PA883393LC01HU
Uniprot NO.	Q9UHD9
Target Names	UBQLN2
Alternative Names	UBQLN2 antibody; N4BP4 antibody; PLIC2 antibody; HRIHFB2157Ubiquilin-2 antibody; Chap1 antibody; DSK2 homolog antibody; Protein linking IAP with cytoskeleton 2 antibody; PLIC-2 antibody; hPLIC-2 antibody; Ubiquitin-like product Chap1/Dsk2 antibody
Raised in	Rabbit
Species Reactivity	Human
Immunogen	Recombinant Human Ubiquilin-2 protein (473-542AA)
Immunogen Species	Homo sapiens (Human)
Conjugate	FITC
Isotype	IgG
Purification Method	>95%, Protein G purified
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Image	 UBQLN2 Antibody, FITC conjugated