




FOSB Antibody, Biotin conjugated

Product Code	CSB-PA008791LD01HU
Uniprot NO.	P53539
Target Names	FOSB
Alternative Names	Activator protein 1 antibody; AP 1 antibody; DKFZp686C0818 antibody; FBJ murine osteosarcoma viral oncogene homolog B antibody; fosB antibody; FOSB_HUMAN antibody; G0/G1 switch regulatory protein 3 antibody; GOS3 antibody; GOS 3 antibody; GOS B antibody; GOS3 antibody; GOSB antibody; MGC42291 antibody; Oncogene FOS B antibody; Oncogene FOSB antibody; Protein fos B antibody; Protein fosB antibody
Raised in	Rabbit
Species Reactivity	Human
Immunogen	Recombinant Human Protein fosB protein (1-41AA)
Immunogen Species	Homo sapiens (Human)
Conjugate	Biotin
Isotype	IgG
Purification Method	>95%, Protein G purified
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Tested Applications	ELISA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Lead Time	Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.
Image	 FOSB Antibody, Biotin conjugated