



AKT1 Antibody

Product Code	CSB-PA612580XA01OFG
Uniprot NO.	Q0JKV1
Target Names	AKT1
Alternative Names	AKT1 antibody; Os01g0648000 antibody; LOC_Os01g45990 antibody; Potassium channel AKT1 antibody; OsAKT1 antibody
Raised in	Rabbit
Species Reactivity	Oryza sativa subsp. japonica (Rice)
Immunogen	Recombinant Oryza sativa subsp. japonica (Rice) AKT1 protein
Conjugate	Non-conjugated
Isotype	IgG
Purification Method	Antigen Affinity Purified
Concentration	It differs from different batches. Please contact us to confirm it.
Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Tested Applications	ELISA, WB (ensure identification of antigen)
Lead Time	Made-to-order (14-16 weeks)
Value-added Deliverables	200ug * antigen (positive control); 1ml * Pre-immune serum (negative control);
Quality Guarantee	Antibody purity can be guaranteed above 90% by SDS-PAGE detection; ELISA titer can be guaranteed 1: 64,000; WB validation with antigen can be guaranteed positive;