




Recombinant Mouse Insulin-like growth factor II (Igf2)

Product Code	CSB-EP011088MO
Purity	Greater than 85% as determined by SDS-PAGE.
Target Names	Igf2
Uniprot NO.	P09535
Research Area	Signal Transduction
Species	Mus musculus(Mouse)
Source	E.coli
Expression Region	25-91aa
Target Protein Sequence	AYGPGETLCGGELVDTLQFVCSDRGFYFSRPSSRANRRSRGIVEECCFRSCDL ALLETYCATPAKSE
Mol. Weight	13.4 kDa
Protein Length	Full Length of Mature Protein
Tag Info	N-terminal 6xHis-tagged
Buffer	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Lead Time	3-7 business days
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.