



# CLTC Antibody, Biotin conjugated

|                            |                                                                                                                                                                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA005593LD01HU                                                                                                                                                                                                                            |
| <b>Uniprot NO.</b>         | Q00610                                                                                                                                                                                                                                        |
| <b>Target Names</b>        | CLTC                                                                                                                                                                                                                                          |
| <b>Alternative Names</b>   | Clathrin heavy chain 1 (Clathrin heavy chain on chromosome 17) (CLH-17), CLTC, CLH17 CLTCL2 KIAA0034                                                                                                                                          |
| <b>Raised in</b>           | Rabbit                                                                                                                                                                                                                                        |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                                                                         |
| <b>Immunogen</b>           | Recombinant Human Clathrin heavy chain 1 protein (851-990aa)                                                                                                                                                                                  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                          |
| <b>Conjugate</b>           | Biotin                                                                                                                                                                                                                                        |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                           |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                                                                                      |
| <b>Concentration</b>       | It differs from different batches. Please contact us to confirm it.                                                                                                                                                                           |
| <b>Buffer</b>              | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                                                                                              |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                                         |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                 |
| <b>Lead Time</b>           | Basically, we can dispatch the products out in 1-3 working days after receiving your orders. Delivery time maybe differs from different purchasing way or location, please kindly consult your local distributors for specific delivery time. |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                                                                   |
| <b>Image</b>               |                                                                                                                                                                                                                                               |

CLTC Antibody, Biotin conjugated

