




Recombinant Human Fibroblast growth factor receptor 4 (FGFR4), partial

Product Code	CSB-MP008648HU1
Abbreviation	Recombinant Human FGFR4 protein, partial
Purity	Greater than 90% as determined by SDS-PAGE.
Target Names	FGFR4
Uniprot NO.	P22455
Research Area	Cancer
Alternative Names	FGFR-4;CD334;FGFR4;JTK2, TKF
Species	Homo sapiens (Human)
Source	Mammalian cell
Expression Region	22-369aa
Target Protein Sequence	LEASEEVELEPCLAPSLEQQEQELTVALGQPVRLCCGRAERGHHWYKEGSRL APAGRVRGWRGRLEIASFLPEDAGRYLCLARGSMIVLQNLTLITGDSLTSND DEDPKSHRDPNHRHSYPQQAPYWTHPQRMEKKLHAVPAGNTVKFRCPAAG NPTPTIRWLKDGQAFHGENRIGGIRLRHQHWSLVMSVVPSPDRGTYTCLVEN AVGSIRYNYLLDVLERSPHRPILQAGLPANTTAVVGSDEVLLCKVYSDAQPHI QWLKHIVINGSSFGADGFPYVQVLKTADINSSEVEVLRLNVSAEDAGEYTC LAGNSIGLSYQSAWLTVLPEEDPTWTAAPEARYTD
Mol. Weight	67.4 kDa
Protein Length	Partial
Tag Info	C-terminal hFc1-tagged
Form	Liquid or Lyophilized powder
Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5% -50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20/-80. The shelf life of lyophilized form is 12 months at -20/-80.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Lead Time	Basically, we can dispatch the products out in 3-7 working days after receiving your orders. Delivery time may differ from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.
