



Recombinant Macaca mulatta (Rhesus macaque) Interleukin-6 (IL6)

Product Code	011664MOW
Purity	>85% (SDS-PAGE)
Target Names	IL6
Uniprot NO.	P51494
Species	Macaca mulatta (Rhesus macaque)
Source	Yeast
Expression Region	28-212
Target Protein Sequence	APV LPGEDSKNVA APHSQPLTSS ERIDKHIRYI LDGISALRKE TCNRSNMCES SKEALAENNL NLPKMAEKDG CFQSGFNEDT CLVKIITGLL EFEVYLEYLQ NREFESSEEQA RAVQMSTKVL IQFLQKKAKN LDAITTPEPT TNASLLTKLQ AQNQWLQDMT THLILRSFKE FLQSNLRALR QM
Protein Length	Full Length of Mature Protein
Tag Info	N-terminal His-tagged/Tag-Free
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.