



Recombinant Pongo abelii OX-2 membrane glycoprotein (CD200), partial

Product Code	722381PYX
Abbreviation	CD200
Purity	>85% (SDS-PAGE)
Target Names	CD200
Uniprot NO.	Q5RAL8
Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Source	Yeast
Protein Length	Partial
Tag Info	N-terminal His-tagged/Tag-Free
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.