



Recombinant Human Interleukin-17 receptor E (IL17RE), partial

Product Code	839874HU
Abbreviation	IL17RE
Purity	>85% (SDS-PAGE)
Target Names	IL17RE
Uniprot NO.	Q8NFR9
Species	Homo sapiens (Human)
Source	Yeast
Protein Length	Partial
Tag Info	N-terminal His-tagged/Tag-Free
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.